# Development of Falling Accident Response System for the Elderly people

### TIMI-RUM

Project Background

**INDEX** 

Project Progress

Results of Project Execution

# 1. Project Background

# NECESSITY OF PROJECT

🤽 연합뉴스 | 2020.02.03. | 네이버뉴스

65세이상 겨울철 **낙상사고** 4년새 52%↑..."주로 집서 일상생활중"

노인 낙상 장소별로는



M 전남매일 2019.04.24.

<sub>질본 "65</sub> 노년기 웰다잉(Well-Dying) 불청객 '**낙상사고**'


노년기 웰다잉(Well-Dying) 불청객 '낙상사고' 김영돈 광주시 소방안전본부 구조구 급과장 1970년에 62.... 버킷리스트 작성, 건강 체크, 고독사 예방, 자원봉사 참여 ...

> The rate of elderly people's fall accidents has increased.

> > However, elderly people living alone often suffer from absence of appropriate response system if accidents occur.

# 2. Project Progress





# IDEA MATERIALIZATION



### **IMAGE RECOGNITION**

Image recognition using YOLO v3

#### SOUND DETECTION

Sound detection using python

### **DJANGO & REACT**


Able to check the result on web page

### STT & TTS


Inquiring to user about the seriousness when falling is detected



# SYSTEM OVERVIEW DIAGRAM



**Falling sound** 



# IMPLEMENTATION PLAN

1

IMAGE RECOGNITION SOUND DETECTING


JUDGMENT OF FALLING ACCIDENT SEND
JUDGMENT
TO REACT

IMPLEMENTATION
WEB WITH REACT

5

Labeling

Labeling two posture (general, falling) using YOLO-MARK








Falling

# IMPLEMENTATION PLAN

### **Training**


Training dataset using YOLO v3

Dataset is composed to 20,000 labeled photo



Loss Graph

# USAGE OF WEBOS COMPONENT



**REACT.JS** 

For web implementation




**CAMERA** 

**MODULE** 

For image recognition & sound detecting

# 3. Results of Project Execution

# TRAINING RESULT





### **TIMI-RUM**

#### Falling Accident Response System for the Elderly people

As the number of single-person elderly households increases these days, there are many cases in which the situation worsens as appropriate measures are not taken in the event of falls.

The project is carried out with the aim of developing a response system that can detect and report in the event of a fall accident.

GO TO SHOW RESULT





### **TIMI-RUM**







### **Level 0 – Nothing Happened**

The state that falling is not detected yet

# WEB PAGE - Result / Accident Occur

#### TIMI-RUM





## Level 1 – Falling Detected

Falling is detected both web cam screen and sound occurrence at the same time

System ask whether help is needed

# **WEB PAGE** - Result / Emergency



### Level 2 - Help is needed

If the person say yes, then in web page, design is changed to red, and alarm is played

# THANK YOU!